


## HOW EFFICIENT ARE YOUR BIO-CERAMIC TIPS?

WASTE LESS, SEAL MORE WITH AVALON BIOMED™



## HIGHER CEMENT CONTENT

COMPETITOR ANALYSIS OF CEMENT CONTENT<sup>2</sup>




Scan here to read journal abstract



*NeoSEALER EZFlo's flowability marries perfectly with the minimal shapes created by SlimShapers. It's easily extruded through 29ga NaviTips and provides an excellent 3D seal.*

**Kirk Coury, DDS, MS**



*The versatility of NeoSEALER Flo allows for a consistent warm obturation without caking or drying out. It's my go to sealer for all cases!*

**Ryan Walsh, DDS, MS**



*After the canal is retro-prepped, I fill the canal with NeoSEALER, then place a "lid" with NeoPUTTY. It works awesome, has great radiographic appearance, and has given me tremendous success over the years.*

**Jacob Butler, DMD, MS**



*NeoSEALER EZFlo is my go-to bioceramic sealer. The NaviTip ensures precise placement, even in challenging canals, while the sealer flows smoothly and maintains integrity when heated. It provides excellent control and navigates lateral canals seamlessly, delivering optimal results with ease and predictability.*

**Erin Duarte, DDS**



**THE ONLY U.S. BASED COMPANY MANUFACTURING PREMIXED BIOACTIVE BIO-CERAMICS UNDER A U.S. PATENT**

Patent No. 12,171,898

# NeoSEALER<sup>®</sup> EZFlo



DESIGNED FOR MINIMALLY INVASIVE ENDODONTICS (MIE)



Consistent flow through a 29ga NaviTip™ 1\*

Rotating wing design allows for tip to be used in any orientation

**FORMULATED WITH MORE BIOACTIVE CEMENT THAN OTHER BIOCERAMIC SEALERS<sup>2</sup>**

Formulated with MORE bioactive cement than competitors for improved setting and solubility.

**BETTER FLOW AT A BETTER PRICE**

~30% less expensive compared to leading competitor.

**VERSATILE**

One sealer to streamline your practice. Suitable for both warm and cold obturation techniques. NeoSEALER EZFlo will NOT change consistency with warm obturation techniques.

<sup>1</sup> Products listed are not registered trademarks of Avalon Biomed  
<sup>2</sup> Data on file

# NeoSEALER<sup>®</sup> Flo

THE MOST VERSATILE AND ECONOMICAL HIGH CEMENT CONTENT BIOCERAMIC SEALER



**Flex Flo Tip™**  
Operate with more precision and efficiency with 81% less product waste with our Flex Flo Tip minimal-waste design.

**30 min**  
working time

**9.5 hrs**  
setting time\*\*

\*\*Sets in vivo in the presence of moisture

**HIGHER BIOACTIVE CEMENT THAN OTHER BIOCERAMIC SEALERS<sup>2</sup>**

NeoSEALER Flo is thoughtfully designed to promote setting and maximize bioactivity.

**81% LESS PRODUCT WASTE WHEN USING FLEX FLO TIP**

Our unique Flex Flo Tip™ design combines a flexible cannula with an integrated stopper attached to a minimal waste barrel for highly accurate placement.

**RESIN-FREE**

NeoSEALER Flo is resin-free for maximum bioactivity and biocompatibility.

NeoSEALER Flo releases calcium and hydroxide ions from the MTA, promoting hydroxyapatite (HA) formation on the MTA surfaces to enhance sealing and support healing.

Other Available Tip Option:



Basic Flex Tip

# NeoPUTTY<sup>®</sup>

PREMIXED ZERO-WASTE BIOACTIVE BIOCERAMIC



Positive displacement plunger tip for zero waste.

**5 hrs**  
setting time\*

\*Sets in vivo in the presence of moisture

**PERFECT CONSISTENCY, NO MIXING, NO WASTE**

NeoPUTTY's uniform, low-tack consistency from beginning to end, with a 3-year shelf life and no dry out between uses, ensures that no product goes to waste.

**UNMATCHED RADIOCAPACITY**

Highest in pre-mixed bioceramic putty class at 8.1 mm Al equivalent.

**IMMEDIATELY WASH-OUT RESISTANT**

Stays where you want it so you can efficiently complete the procedure ensuring bioactive benefits.

# NeoMTA<sup>®</sup> 2

ONE PRODUCT FOR 13 OF YOUR MOST IMPORTANT PROCEDURES



Eliminate guesswork and reduce product waste with this 0.1 g measurement scoop.

**PUTTY**  
**14 min**  
working & setting time

**SEALER**  
**20 min - 1 hr**  
working & setting time

**ONE PRODUCT FOR EVERY PROCEDURE**

NeoMTA 2's unique proprietary gel allows you to mix to a firm putty or a thinner sealer.

**HIGHER RADIOCAPACITY**

NeoMTA 2 - 6.5 mm Al equivalent at putty consistency. Competitors - 5.0 mm Al equivalent or less.

**IMMEDIATELY WASH-OUT RESISTANT**

When mixed to a putty, NeoMTA 2 stays where you want it so you can efficiently complete the procedure ensuring bioactive benefits.

\*Products listed are not registered trademarks of Avalon Biomed.