



UNIVERSAL CHAIRSIDE SPACE MAINTAINER

Redesigned with YOU in mind



Ask about our Buyback program!

UNIVERSAL DESIGN

= HALF THE INVENTORY

$\frac{1}{2}$

Band designed for both maxillary and mandibular tooth anatomy, reduces inventory by one half



Available in a full range of sizes: 23-42 with $\frac{1}{2}$ sizes



Easier re-ordering process



304L stainless-steel band that is strong and tempered to be more malleable

SHORTER TUBE DESIGN

SAVES VALUABLE CHAIR TIME

No trimming or cleaning of metal debris required, saving you valuable chair time

Tube and loop diameter constructed with a slightly snug fit to minimize dislodgement and breakage

NuSmile

vs.

Other Brands




4 mm

5 mm

STRONG

LASER WELD

Strong laser welded joint for enhanced durability



PREFORMED VS. LAB SPACE MAINTAINERS

Advantages Vs. Lab Fabricated Space Maintainers

1/3rd THE COST
of lab fabricated
space maintainers


**SINGLE
APPOINTMENT**
Clinicians can
fabricate a space
maintainer in just
minutes.
*No Second
Appointment Needed!*

NO IMPRESSION
scanning or model
pouring required

NO REMAKES
Eliminates the
possibility of
laboratory
fabrication errors

WHY MAKE


THE SWITCH



Half the inventory with a
universal band design for
both maxillary and
mandibular tooth anatomy



Higher quality
at a lower cost



Shorter tube design,
saves valuable chair time

SINGLE WIDTH WIRES/LOOPS

NEW


Loops now include 1mm markings for easy and precise trimming

Wires/loops are available in only one width, reducing inventory

Optimal curvature and height to reduce slip and tissue impingement on primary or permanent cuspids

Curved loop can be placed distally in a gingival direction to guide an erupting permanent molar



Distal shoe is designed with slight curvature following proximal molar anatomy, reducing plaque traps



NuSmile

vs.

Other Brands



BUT WAIT!

THERE'S MORE...

Higher quality and less cost than other brands

50% less inventory with universal band design in a full range of sizes 23 to 42 with 1/2 sizes

Shorter tube to reduce time spent trimming and adjusting



Three standard loop options to cover all clinical cases

Tube and loop diameter are constructed with a slightly snug fit to minimize dislodgement and breakage


Universal band is easily adaptable to both maxillary and mandibular tooth anatomy

UNIVERSAL BAND, SHORT TUBE DESIGN

REDUCES INVENTORY & SAVES TIME


NuSmile Space Maintainers provide a convenient solution to minimize space drift after the premature loss of primary molars at a third of the cost of a lab fabricated space maintainer

PLAIN LOOP




- A versatile loop option to suit a variety of clinical scenarios
- 316L loop wires are strong and can be easily adjusted
- 1mm markings for easy and precise trimming

CURVED LOOP



- Optimal curvature and height to reduce slip and tissue impingement on primary or permanent cusps
- Can be placed distally in a gingival direction to guide an erupting permanent molar

DISTAL SHOE




- Used sub-gingivally to guide an erupting permanent molar
- Designed with slight curvature following proximal molar anatomy, reducing plaque traps

NEW

KITS


STARTER KIT

Contains 39 Bands, 18 Plain Loop, 18 Curved Loop, 3 Distal Shoe



PROFESSIONAL KIT

Contains 78 Bands, 36 Plain Loop, 36 Curved Loop, 6 Distal Shoe



REFILL OPTIONS SELECT FROM

1 Pack Band | Plain Loop 10 Pack | Curved Loop 10 Pack | Distal Shoe 10 Pack

QUOTES AND CASES

When my lab increased the price of making a band and loop, I decided to look into the NuSmile Chairside Space Maintainer kit. The price point was more than competitive. I never thought a distal shoe would be so easy to place.



Stephen E Greenleaf, DMD

The NuSmile Space Maintainer with its universal band size and shorter tube makes for less stuff, less fuss, and less chair time. That's a win, win, win for my patients, my team, and me!



Ben Curtis, DDS

I used a competitor's chairside space maintainers for over 20 years, but found that I needed to trim about 30-40% of the tubes in order to make the spacer fit the space. With the NuSmile space maintainers I have not had to trim a single one!



William "Bill" Waggoner, DDS, MS

The NuSmile Universal Chairside Space Maintainer simplifies our practice with its universal band design, cutting inventory needs in half while fitting both maxillary and mandibular molars perfectly. The laser-marked loops and shorter tubes make placement a breeze, ensuring reliable results every time.



David Harney, DDS, MS

